

erccdashboard Package Vignette

Sarah A. Munro

December 15, 2013

This vignette describes the use of the `erccdashboard` R package to analyze External RNA Control Consortium (ERCC) spike-in control ratio mixtures in gene expression experiments. Two types of data from the SEQC/MAQC III project were analyzed.

1. Rat toxicogenomics treatment and control samples for different drug treatments
2. Human reference RNA samples from the MAQC I project, Universal Human Reference RNA (UHRR) and Human Brain Reference RNA (HBRR)

1 Rat Toxicogenomics Example: MET treatment

1.1 Load data and define input parameters

Load the `testerccdashboard` package.

```
> library( "erccdashboard" )
```

Load the Rat Toxicogenomics Data set.

```
> load(file = system.file("data/SEQC.RatTox.Example.RData",  
                           package = "erccdashboard"))
```

The R workspace should now contain 5 count tables and for each count table a corresponding total reads vector. Take a look at the data for the MET experiment.

```
> head(COH.RatTox.ILM.MET.CTL.countTable)
  Feature MET_1 MET_2 MET_3 CTL_1 CTL_2 CTL_3
1 ERCC-00002 16629 18798 26568 36600 45436 25163
2 ERCC-00003  1347  1565  1983  3048  3447  2195
3 ERCC-00004  4569  5570  6755  1240  1484   902
4 ERCC-00009   811   869  1123   909  1073   537
5 ERCC-00013    3    1    2    1    5    1
6 ERCC-00019   24   32   43    5   13    4
> COH.RatTox.ILM.MET.CTL.totalReads
[1] 41423502 46016148 44320280 38400362 47511484 33910098
```

The first column of the count table, `Feature`, contains unique names for all the transcripts that were quantified in this experiment. The remaining columns represent replicates of the pair of samples, in this count table the control sample is labeled CTL and the treatment sample is labeled MET. An underscore is included to separate the sample names from the replicate numbers during analysis. This naming convention `Sample_Rep` is needed for the columns of any input count table.

The total reads vectors will be used for library size normalization of the count tables. Total reads can either represent the total number of reads in FASTQ files or total mapped reads. In the examples provided with this package FASTQ file total reads are used.

For our analysis of the MET-CTL experiment start by assigning the MET-CTL data to the input data variables `countTable` and `totalReads`.

```
> countTable <- COH.RatTox.ILM.MET.CTL.countTable
> totalReads <- COH.RatTox.ILM.MET.CTL.totalReads
```

In addition to `countTable` and `totalReads`, there are 7 additional variables that must be defined by the user. First the filename prefix for results files, `filenameRoot`, needs to be defined. Here we choose to use the lab abbreviation COH and the platform abbreviation ILM as our identifiers, but this is flexible for the user.

```
> filenameRoot = "COH.ILM"
```

Next, 6 parameters associated with the ERCC control ratio mixtures need to be defined, `sample1Name`, `sample2Name`, `erccMixes`, `ERCCdilution`, `spikeVol`, and `totalRNAmass`.

The sample spiked with ERCC Mix 1 is `sample1Name` and the sample spiked with ERCC Mix 2 is `sample2Name`. In this experiment `sample1Name` = MET and `sample2Name` = CTL. For a more robust experimental design the reverse spike-in design could be created using additional replicates of the treatment and control samples. ERCC Mix 2 would be spiked into MET samples and ERCC Mix 1 would be spiked into CTL control replicates.

`ERCCdilution` is the dilution factor of the pure Ambion ERCC mixes prior to spiking into total RNA samples. Here a 1/100 dilution was made from the Ambion ERCC mixes according to the protocol. The amount of diluted ERCC mix spiked into the total RNA sample is `spikeVol` (units are μL). The mass of total RNA spiked with the diluted ERCC mix is `totalRNAmass` (units are μg).

```
> sample1Name = "MET"
> sample2Name = "CTL"
> ERCCMixes = "Ambion4plexPair"
> ERCCdilution = 1/100
> spikeVol = 1
> totalRNAmass = 0.500
```

The final required input parameter, `choseFDR`, is the False Discovery Rate (FDR) for differential expression testing. A typical choice would be 0.05 (5% FDR), for the rat data sets a more liberal FDR was used, `choseFDR` = 0.1.

```
> choseFDR = 0.1
```

In addition to the required input variables the user can also choose whether to print the results directly to a PDF file (the default is TRUE) with the variable `printPDF`.

1.2 Initialize the `expDat` list for analysis

The `expDat` list is created with the `initDat` function:

```
> expDat <- initDat(countTable, totalReads, filenameRoot, sample1Name,
                    sample2Name, ERCCMixes, ERCCdilution, spikeVol, totalRNAmass,
                    choseFDR)
```

Filename root is: COH.ILM.MET.CTL

Library sizes:

41.4235 46.01615 44.32028 38.40036 47.51148 33.9101

Using total sequencing reads,

mean library size factor = 41.93031

Look at the structure of `expDat`

```
> summary(expDat)
```

	Length	Class	Mode
sampleInfo	16	-none-	list
totalReads	6	-none-	numeric
Transcripts	7	data.frame	list
designMat	3	data.frame	list
sampleNames	2	-none-	character
idCols	6	data.frame	list
totalReads	6	-none-	numeric
expressDat	7	data.frame	list
libeSize	6	-none-	numeric
ERCCxlabelIndiv	1	-none-	expression
ERCCxlabelAve	1	-none-	expression
spikeFraction	1	-none-	numeric
mnLibeFactor	1	-none-	numeric
sampleLibeSums	6	-none-	numeric

The expDat list will be passed to the erccdashboard functions for analysis of technical performance.

1.3 Estimate the mRNA fraction difference, r_m for the pair of samples

Estimate r_m for the sample pair using a negative binomial glm. The r_m results will be added to the expDat structure and are necessary for the remaining analysis.

```
> expDat <- est_r_m(expDat)
Check for sample mRNA fraction differences(r_m)...

log.offset
17.53936 17.6445 17.60695 17.46358 17.67648 17.33922

Number of ERCC Controls Used in r_m estimate
63

Outlier ERCCs for GLM r_m Estimate:
None

GLM log(r_m) estimate:
-0.0472291

GLM log(r_m) estimate standard deviation:
0.02061546

GLM r_m estimate:
0.9538688

GLM r_m upper limit
0.9599503

GLM r_m lower limit
0.9478259
```

An r_m of 1 indicates that the two sample types under comparison have similar mRNA fractions of total RNA. The r_m estimate is used to adjusted the expected ERCC mixture ratios in this analysis and may indicate a need for a different sample normalization approach.

1.4 Test for differential expression

Test for differential expression with the `geneExprTest` function. This function wraps the QuasiSeq differential expression testing package. If a correctly formatted csv file is provided with the necessary DE test results, then `geneExprTest` will bypass DE testing (with reduced runtime). The function will look for a csv file with the name "filenameRoot.quasiSeq.res.csv" and the first 3 column headers of the file must be "Feature", "pvals", and "qvals".

```
> expDat <- geneExprTest(expDat)
```

1.5 Diagnostic Performance: ROC curves and AUC statistics

Generate ROC curves for the differential ratios and the corresponding Area Under the Curve (AUC) statistics.

```
> expDat = erccROC(expDat)
Area Under the Curve (AUC) Results:
  Ratio  AUC Measured Spiked
1   4:1 1.000      16      23
2 1:1.5 0.950      16      23
3   1:2 0.967      16      23
```

1.6 Diagnostic Performance: Limit of Detection of Ratios (LODR)

Find LODR estimates using the ERCC data p-values.

```
> expDat = estLODR(expDat, kind = "ERCC", prob=0.9)
  Feature  Feature  MnCnt  Pval  LogPval
1      1 1 ERCC-00002 28199.000000 2.283076e-03 -6.0822315
2      2 2 ERCC-00003 2264.166667 4.464509e-04 -7.7141812
3      3 3 ERCC-00004 3420.000000 5.301011e-06 -12.1476129
4      4 4 ERCC-00009  887.000000 8.860051e-01 -0.1210326
5      5 5 ERCC-00013   2.166667 6.627015e-01 -0.4114307
6      6 6 ERCC-00019  20.166667 1.703913e-03 -6.3748279
  F.stat Fold  Den.df
1 20.39539387 0.5 3.531005
2 35.03013540 0.5 3.531005
3 128.67218955 4.0 3.531005
4  0.02199068 1.0 3.531005
5  0.20601552 0.5 3.531005
6 22.57265968 4.0 3.531005
  Ratio LODR Estimate 90% CI Lower Bound 90% CI Upper Bound
1   4:1           26           19           32
3 1:1.5           Inf          <NA>          <NA>
4   1:2          250          120          350
```

One can also obtain LODR estimates using p-values simulated from endogenous transcripts

```
> expDat = estLODR(expDat, kind = "Sim", prob = 0.9)
```

	Feature	Feature	MnCnt	Pval	LogPval
1	Sim0.5Fold1	Sim0.5Fold1	2.333333	0.01660461	-4.0980748
2	Sim0.5Fold2	Sim0.5Fold2	2.166667	0.31560417	-1.1532665
3	Sim0.5Fold3	Sim0.5Fold3	2.500000	0.57873930	-0.5469032
4	Sim0.5Fold4	Sim0.5Fold4	4.000000	0.59722891	-0.5154548
5	Sim0.5Fold5	Sim0.5Fold5	3.833333	0.80939536	-0.2114678
6	Sim0.5Fold6	Sim0.5Fold6	10.666667	0.03936595	-3.2348539

	F.stat	Fold	Den.df
1	9.38122816	0.5	3.561428
2	1.15497020	0.5	3.561428
3	0.33647681	0.5	3.561428
4	0.30411297	0.5	3.561428
5	0.06239896	0.5	3.561428
6	6.18334070	0.5	3.561428

	Ratio	LQDR	Estimate	90% CI Lower Bound	90% CI Upper Bound
1	4:1		35	21	41
3	1:1.5		Inf	<NA>	<NA>
4	1:2		Inf	<NA>	<NA>

1.7 Use dynRangePlot function to evaluate dynamic range of the data

This function will add a plot to `expDat$Figures` of the signal vs. abundance of the spiked ERCC controls.


```
> expDat <- dynRangePlot(expDat, errorBars = T)
Number of ERCCs in Mix 1 dyn range: 63

Number of ERCCs in Mix 2 dyn range: 63
These ERCCs were not included in the signal-abundance plot,
because not enough non-zero replicate measurements of these
controls were obtained for both samples:
```

ERCC-00058
ERCC-00067
ERCC-00077
ERCC-00168
ERCC-00028
ERCC-00033
ERCC-00040
ERCC-00109
ERCC-00154
ERCC-00158

View the plot that is now stored in `expDat` with the command

```
> expDat$Figures$plotdynRange
```



This figure shows that in this experiment the expected signal-abundance relationship spans a 2^{15} dynamic range. To capture the full 2^{20} dynamic range design of the control mixtures additional sequencing depth may be needed.

1.8 Use LODR estimates to Annotate Signal-Abundance and Ratio-Abundance Plots

Get LODR annotations for adding to plots and then annotate the signal-abundance and ratio-abundance plots using LODR estimate information

```
> expDat <- annotLODR(expDat)
```

	Fold	Ratio	Count	Log2Count_normalized	Log2Conc
1	4:1	4:1	26	-0.689482	-1.106341
2	1:1	1:1	NA	NA	NA
3	1:1.5	1:1.5	Inf	Inf	Inf
4	1:2	1:2	250	2.575863	2.232469

LODR estimates are available to code ratio-abundance plot

These ERCCs were not included in the ratio-abundance plot,
because not enough non-zero replicate measurements of these
controls were obtained for both samples:

ERCC-00028

ERCC-00033

ERCC-00040

ERCC-00058

ERCC-00067

ERCC-00077

ERCC-00109

ERCC-00154

ERCC-00158

ERCC-00168

Global Ratio SD for this sample pair is: 0.7785041

Estimate Std. Error

Minimum SD Estimate 0.4787852 0.03264151


Maximum SD Estimate 1.6195817 0.20936506

Lambda 0.4123437 0.05767326


1.9 Viewing Diagnostic Plots

All dashboard plots are stored in the `expDat$Figures` list. You can call any figure for viewing directly and you can also save the figures to a pdf file. The six plots presented in the `erccd` dashboard publication can be generated with the following commands.

```
> expDat$Figures$ROCplot
```





```
> expDat$Figures$plotLODR.ERCC
```




Ratio	LODR Estimate	90% CI Lower Bound	90% CI Upper Bound
4:1	26	19	32
1:1.5	Inf	NA	NA
1:2	250	120	350


```
> expDat$Figures$dispPlot
```




```
> expDat$Figures$plotdynRangeAnnot
```



```
> expDat$Figures$plotERCCeffects
```



```
> expDat$Figures$plotRatioAnnot
```



The function `savePlots` will save selected figures to a pdf file. The default is the 6 manuscript figures to a single page (`plotsPerPg = "manuscript"`). If `plotsPerPg = "single"` then each plot is placed on an individual page in one pdf file. If `plotlist` is not defined (`plotlist = NULL`) then all plots in `expDat$Figures` are printed to the file.

```
> savePlots(expDat)
```

1.10 Output Results for Comparisons Across Experiments or Between Laboratories

If you wish, save your results to an Rdata file that can be reused for comparisons across experiments or between laboratories.

```
> saveResults(expDat)
```

2 SEQC Reference RNA Examples: UHRR vs. HBRR

2.1 Load data and define input parameters

```
> load(file = system.file("data/SEQC.Main.Example.RData",
                           package = "erccdashboard"))
> countTable <- Lab5.ILM.UHRR.HBRR.countTable
> totalReads <- Lab5.ILM.UHRR.HBRR.totalReads
> filenameRoot = "Lab5"
> sample1Name = "UHRR"
> sample2Name = "HBRR"
> ERCCMixes = "Ambion4plexPair"
> ERCCdilution = 1
> spikeVol = 50
> totalRNAmass = 2.5*10^(3)
> choseFDR = 0.1
> expDat <- initDat(countTable, totalReads, filenameRoot, sample1Name,
                    sample2Name, ERCCMixes, ERCCdilution, spikeVol, totalRNAmass,
                    choseFDR)
```

Filename root is: Lab5.UHRR.HBRR

Library sizes:

138.7869 256.0065 199.4683 431.9338 247.9856 219.3833 251.2658 257.5082

Using total sequencing reads,

mean library size factor = 250.2923

```
> expDat <- est_r_m(expDat)
```

Check for sample mRNA fraction differences(r_m)...

log.offset

18.74845 19.36071 19.11117 19.88378 19.32888 19.20633 19.34202 19.36656

Number of ERCC Controls Used in r_m estimate

71

```

Outlier ERCCs for GLM r_m Estimate:
ERCC-00137 ERCC-00085 ERCC-00054 ERCC-00019

```

```

GLM log(r_m) estimate:
0.2335326

```

```

GLM log(r_m) estimate standard deviation:
0.002941665

```

```

GLM r_m estimate:
1.263054

```

```

GLM r_m upper limit
1.264133

```

```

GLM r_m lower limit
1.261975

```

```

> expDat <- geneExprTest(expDat)

```

2.2 Diagnostic Performance: ROC curves and AUC statistics

Generate ROC curves for the differential ratios and the corresponding Area Under the Curve (AUC) statistics.

```

> expDat = erccROC(expDat)
Area Under the Curve (AUC) Results:
  Ratio  AUC Measured Spiked
1   4:1 1.000      17      23
2  1:1.5 0.967      17      23
3   1:2 0.991      19      23

```

2.3 Diagnostic Performance: Limit of Detection of Ratios (LODR)

Find LODR estimates using the ERCC data p-values.

```

> expDat = estLODR(expDat, kind = "ERCC", prob=0.9)
  Feature  Feature  MnCnt  Pval  LogPval
1      1 ERCC-00002 1162406.125 2.803132e-22 -49.6261347
2      2 ERCC-00003  95757.625 2.449837e-25 -56.6686058
3      3 ERCC-00004 128395.250 1.194356e-32 -73.5051157
4      4 ERCC-00009  36223.125 5.248996e-01  -0.6445483
5      5 ERCC-00013   97.625 4.840217e-06 -12.2385511
6      6 ERCC-00014   67.125 1.948995e-04  -8.5430264
  F.stat Fold  Den.df
1  874.1474841  0.5 21.31283
2 1481.8037744  0.5 21.31283
3 5148.2427909  4.0 21.31283
4   0.4148318  1.0 21.31283
5  32.1734341  0.5 21.31283

```

```

6 18.4971302 0.5 21.31283
Ratio LODR Estimate 90% CI Lower Bound 90% CI Upper Bound
1 4:1 <27 <27 32
3 1:1.5 61 <25 77
4 1:2 26 <19 30

```

One can also obtain LODR estimates using p-values simulated from endogenous transcripts

```

> expDat = estLODR(expDat, kind = "Sim", prob = 0.9)
      Feature      Feature      MnCnt      Pval      LogPval
1 Sim0.5Fold1 Sim0.5Fold1 25.250 6.017069e-03 -5.113155
2 Sim0.5Fold2 Sim0.5Fold2 45.625 1.670471e-04 -8.697235
3 Sim0.5Fold3 Sim0.5Fold3 50.375 6.023930e-03 -5.112015
4 Sim0.5Fold4 Sim0.5Fold4 77.875 8.617227e-06 -11.661747
5 Sim0.5Fold5 Sim0.5Fold5 110.250 1.320845e-05 -11.234654
6 Sim0.5Fold6 Sim0.5Fold6 140.125 1.214670e-05 -11.318453
      F.stat Fold      Den.df
1 8.860273 0.5 21.49489
2 18.955606 0.5 21.49489
3 8.857510 0.5 21.49489
4 29.694501 0.5 21.49489
5 27.996365 0.5 21.49489
6 28.325308 0.5 21.49489
Ratio LODR Estimate 90% CI Lower Bound 90% CI Upper Bound
1 4:1 <23 <23 <23
3 1:1.5 30 <28 50
4 1:2 <25 <25 27

```

2.4 Use dynRangePlot function to evaluate dynamic range of the data

This function will add a plot to `expDat$Figures` of the signal vs. abundance of the spiked ERCC controls.

```


> expDat <- dynRangePlot(expDat,errorBars=T)
Number of ERCCs in Mix 1 dyn range: 71

Number of ERCCs in Mix 2 dyn range: 71

```

View the plot that is now stored in `expDat` with the command


```
> expDat$Figures$plotdynRange
```



This figure shows that in this experiment the expected signal-abundance relationship spans the full 2^{20} dynamic range design of the control mixtures.

2.5 Use LODR estimates to Annotate Signal-Abundance and Ratio-Abundance Plots

Get LODR annotations for adding to plots and then annotate the signal-abundance and ratio-abundance plots using LODR estimate information

```
> expDat <- annotLODR(expDat)
```

	Fold	Ratio	Count	Log2Count_normalized	Log2Conc
1	4:1	4:1	27	-3.212583	-5.922734
2	1:1	1:1	NA	NA	NA
3	1:1.5	1:1.5	61	-2.036733	-4.778131
4	1:2	1:2	26	-3.267030	-5.975735


LODR estimates are available to code ratio-abundance plot

Global Ratio SD for this sample pair is: 0.3571987
Estimate Std. Error
Minimum SD Estimate 0.04707465 0.005041491
Maximum SD Estimate 0.51313956 0.038070106
Lambda 0.34093086 0.050686966


2.6 Viewing Diagnostic Plots

All dashboard plots are stored in the `expDat$Figures` list. You can call any figure for viewing directly and you can also save the figures to a pdf file. The six plots presented in the `erccd` dashboard publication can be generated with the following commands.

```
> expDat$Figures$ROCplot
```




```
> expDat$Figures$plotLODR.ERCC
```




Ratio	LODR Estimate	90% CI Lower Bound	90% CI Upper Bound
4:1	<27	<27	32
1:1.5	61	<25	77
1:2	26	<19	30


```
> expDat$Figures$dispPlot
```




```
> expDat$Figures$plotdynRangeAnnot
```



```
> expDat$Figures$plotERCCeffects
```



```
> expDat$Figures$plotRatioAnnot
```



The function `savePlots` will save selected figures to a pdf file. The default is the 6 manuscript figures to a single page (`plotsPerPg = "manuscript"`). If `plotsPerPg = "single"` then each plot is placed on an individual page in one pdf file. If `plotlist` is not defined (`plotlist = NULL`) then all plots in `expDat$Figures` are printed to the file.

```
> savePlots(expDat)
```

2.7 Output Results for Comparisons Across Experiments or Between Laboratories

If you wish, save your results to an Rdata file that can be reused for comparisons across experiments or between laboratories.

```
> saveResults(expDat)
```